



Hemo ID™

DONOR QUICK SCREEN (DQS) PANEL

Cost-effective extended antigen typing for rare donor identification

FACILITATING ACCESS TO RARE DONORS

Every year, international blood donor programs lose rare donors due to sickness, age, or other factors. The Hemo ID™ DQS Panel, developed for use on the MassARRAY® System, simplifies the identification process for new donors by quickly and easily testing thousands of samples.

AN IDEAL LABORATORY WORKFLOW

With a rapid and robust protocol that takes only five and a half hours from DNA to report, you can obtain predicted phenotypes from as few as 12 to as many as 2,000 samples per day.

BLOOD GROUP ANTIGENS INCLUDED IN THE HEMO ID™ DQS PANEL

Rh: ce, Ce, cE, CE, C^w, C^x, V, VS, hrB, hrS, Crawford

Kell: k/ K, Kp^a/ Kp^b, Js^a/ Js^b

Duffy: Fy^a/ Fy^b, Fy^x, GATA

Kidd: Jk^a/ Jk^b

MNS: M/ N, S/ s, U, U^{var} P2, U^{var} NY, U^{neg}

Lutheran: Lu^a/ Lu^b

Dombrock: Do^a/ Do^b, Hy +/-, Jo^a +/-

Landsteiner Wiener: LW^a/ LW^b

Diego: Di^a, Di^b

Colton: Co^a/ Co^b

Scianna: Sc1/ Sc2

Cartwright: Yt^a/ Yt^b

BETA GLOBIN VARIANTS

Hemoglobin C: HbC

Hemoglobin E: HbE

Hemoglobin S: HbS

See details on reverse side

For Research Use Only. Not for use in diagnostic procedures.

Agena Bioscience, Inc.
4755 Eastgate Mall
San Diego, CA 92121
Phone: +1.858.882.2800

Order Desk: +1.858.202.9301
Order Desk Fax: +1.858.202.9220
orderdesk@AgenaBio.com
Web: agenabioscience.com

US: +1.877.4.GENOME
EU: +49.40.899676.0
AP: +61.7.3088.1600
CN: +86.21.6427.0566 51-20071R4.0 1216

Agena
BIOSCIENCE

HOW IT WORKS

Genomic DNA samples are amplified using the provided PCR primers, PCR enzyme, and dNTPs. Residual dNTPs are dephosphorylated with SAP. Following single nucleotide primer extension and termination reactions, the extension products are transferred to a support surface (SpectroCHIP® Array) and loaded into the MassARRAY mass spectrometer for data acquisition. Data analysis software then makes genotype assignments and generates predicted phenotypes.

PERFORMANCE

Genotype and predicted phenotype call rates of > 99% were observed at two independent test sites (n= 768). There was a greater than 99% concordance between the predicted phenotype calls generated by the Hemo ID DQS Panel and those derived from using the Immucor PreciseType™ HEA Test or serology (n=88).

ORDERING INFORMATION

CAT NO	FORMAT	SAMPLES/KIT
17934	10x96	192

REFERENCES

- Gassner C, Meyer S, Frey BM, Vollmert C. Matrix-assisted laser desorption/ionisation, time-of-flight mass spectrometry-based blood group genotyping-- the alternative approach. *Transfus Med Rev.* 2013;27(1):2-9.
- Meyer S, Vollmert C, Trost N, Brönnimann C, Gottschalk J, Andreas B, Frey BM, and Gassner C. High-throughput Kell, Kidd, and Duffy matrix-assisted laser desorption/ionization, time-of-flight mass spectrometry-based blood group genotyping of 4000 donors shows close to full concordance with serotyping and detects new alleles. *Transfusion.* 2014;54(12):3198-207.

HEMO ID™ DQS PANEL GENES AND NUCLEOTIDES

BLOOD GROUP	ANTIGENS	GENE	NUCLEOTIDE
Rh	ce, Ce, cE, CE	RHCE	307C>T, i2(109bp ins), 676G>C
	C ^w	RHCE	122A>G
	C ^x	RHCE	106G>A
	V / VS / hrB	RHCE	48G>C, 733C>G, 1006G>T
	hrS	RHCE	48G>C, 712A>G, 733 C>G
	Crawford	RHCE	48G>C, 697C>G, 733C>G
Kell	k / K	KEL	578C>T
	Kp ^a / Kp ^b	KEL	841C>T
	Js ^a / Js ^b	KEL	1790T>C
Duffy	Fy ^a / Fy ^b	DARC	125A>G
	Fy ^x	DARC	125A>G, 265 C>T
	GATA	DARC	-67T>C, 125A>G
Kidd	Jk ^a / Jk ^b	SLC14A1	838A>G
MNS	M / N	GYP A	59C>T; 72T>G
	S / s	GYP B	143C>T
	U	GYP B	Intron 5 +5G>T
	U ^{var} P2	GYP B	143C>T; Intron5+5G>T
	U ^{var} NY	GYP B	143C>T; 208G>T; 230C>T; 251C>G
	U ^{neg}	GYP B	Deleted GYP B
Lutheran	Lu ^a / Lu ^b	BCAM	230G>A
Dombrock	Do ^a / Do ^b	ART4	793A>G
	Hy+/-	ART4	323G>T
	Jo ^a +/-	ART4	350C>T
Landsteiner Wiener	LW ^a / LW ^b	ICAM4	299A>G
Diego	Dj ^a / Dj ^b	SLC4A1	2561C>T
Colton	Co ^a / Co ^b	AQP1	134C>T
Scianna	Sc1 / Sc2	ERMAP	169G>A
Cartwright	Yt ^a / Yt ^b	ACHE	1057C>A
HEMOGLOBIN TYPE	GENE	NUCLEOTIDE	
Hemoglobin C (HbC)	HBB	19G>A	
Hemoglobin E (HbE)	HBB	79G>A	
Hemoglobin S (HbS)	HBB	20A>T	

vB

For Research Use Only. Not for use in diagnostic procedures.

MassARRAY and SpectroCHIP are registered trademarks of Agena Bioscience, Inc. Hemo ID and Agena Bioscience are trademarks of Agena Bioscience, Inc. All other trademarks are owned by their respective companies.
©2016 Agena Bioscience, Inc. All rights reserved.