

SAMPLE COMPOSITION



Chimerism Testing

For Bone Marrow
Engraftment Monitoring

For Research Use Only.
Not for use in diagnostic procedures.

Agena
BIOSCIENCE

Chimerism Testing

Bone Marrow Engraftment Monitoring

Recipients of allogeneic bone marrow engraftments require clinical monitoring to allow for early diagnosis of such post-transplant adverse events as rejection, graft vs. host disease or malignancy relapse. Chimerism testing is performed on specimens to determine the genetic contribution from the transplant recipient and the donor. Short tandem repeat (STR) analysis has traditionally been the standard of care for this testing.

The Challenge for Labs:

Approximately 20,000 bone marrow engraftments are performed in the US each year². In Europe, there are more than 40,000 annually³. This number is expected to grow over the next 10 years and the prevalence of survivors will grow five-fold through 2030^{2,4}. Increased transplant numbers and longer survival mean greater demand for bone marrow chimerism testing.

Growing Demand Strains Lab Resources

Lab resources are strained by increasing sample volumes that must be processed using an inefficient technology.

Inefficient Testing Methods

STR-based bone marrow engraftment chimerism testing is burdensome and requires time-intensive results interpretation by specialized staff.



Time-consuming results calculations burden laboratory resources. **Is there a way to streamline chimerism testing** to save time and money?

A More Efficient Method

Agena's SNP-based Chimeric ID panel and automated results reporting reduces the time and money required to determine composition of a DNA sample.



Automation Friendly Workflow

Results in an 8-hour shift with minimal hands-on time.
No DNA dilution required.



Accurate Chimerism Determination

Accurately determine percent
DNA contribution.



Automated Results Reporting

Chimerism results are calculated in seconds.
Save time with reporting software provided at
no extra cost.



Save Time & Money

Low-cost reagents and automated results reporting saves
lab resources.



Versatile Application Offering

The MassARRAY system enables testing in multiple
application areas including oncology, PGx and
hereditary genetic testing.





Chimeric ID Panel & Results Reporter

The Chimeric ID panel and results reporter for the MassARRAY system provides simplified chimerism testing. The panel consists of 92 independent SNPs with minor allele frequency between 0.45 and 0.55 across several HapMap populations. The results reporter automatically queries stored reference profiles and calculates the composition of post-transplant follow-up specimens in seconds. Results are provided in a convenient report.

ACCURATE CHIMERISM CALLING

Studies demonstrate accurate determination of major and minor DNA contribution in mixed samples across a range of dilutions. Chimerism calls at all dilution points averaged +/- 1% variance from "truth" as established by X/Y chromosome analysis.



NO DNA DILUTION REQUIRED

STR-based methods warn that high DNA input can affect accuracy. Therefore, many labs tediously dilute extracted DNA. The Chimeric ID panel provides accurate chimerism calling with 5ng – 260ng of DNA input. Skip the DNA dilution and save time!



AUTOMATED RESULTS CALCULATION

Percent contribution from major and minor profiles is calculated in seconds and displayed in an easy-to-interpret sample report. This reduces the burden on specialized staff and frees up critical resources.

Group ID: 878		BME Chimerism ID Software Version: "1.36"	
Sample ID	Chimerism Call		
ID: 434001209 Date: 2018-08-30 15:33:01	Chimeric		
Result: Chimeric/Non-Chimeric		Warnings: NA	
Donor ID: 100_F1	83.5%		
Recipient ID: 100_M3	16.5%		

Group ID: 123		BME Chimerism ID Software Version: "1.36"	
MIXED ID: D17.12574 SCT1 rep Date: 2018-08-30 15:33:01	MIXED ID: D17.12574 SCT1 Date: 2018-08-30 15:33:01	DONOR ID: D17.12570 donor 91.2%	DONOR ID: D17.12570 donor 91.4%
RECIPIENT ID: D17.12573 HOST 8.8%	RECIPIENT ID: D17.12573 HOST 8.6%		
MIXED ID: D17.12574 donor Date: 2018-08-30 15:33:01	MIXED ID: D17.12576 SCT2 Date: 2018-08-30 15:33:01	DONOR ID: D17.12570 donor 100%	DONOR ID: D17.12570 donor 86.4%
RECIPIENT ID: D17.12573 HOST 8.80%	RECIPIENT ID: D17.12573 HOST 13.6%		

HISTORIC RESULTS TRACKING

Monitor results over time. All previous results can be instantly displayed in a convenient report.

INTUITIVE SAMPLE LOG-IN

The intuitive sample log-in interface allows you to enter sample names and quickly identify them as post-transplant follow-up, pre-transplant profile, or control samples. A group ID automatically links follow-up samples to the appropriate profiles.

Agena BIOSCIENCE							
BME Chimerism Report Relation Entry							
Identify sample type and Group Identifier.							
Sample ID	Mixed	Recipient	Donor	NTC	Control	Group ID	MultiDonor
D17.12570 donor	<input type="radio"/>	<input type="radio"/>	<input checked="" type="radio"/>	<input type="radio"/>	<input type="radio"/>		<input type="checkbox"/>
D17.12573 host	<input type="radio"/>	<input checked="" type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>		<input type="checkbox"/>
D17.12574 SCT1	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>		<input type="checkbox"/>
D17.12574 SCT1 rep	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>		<input type="checkbox"/>
D17.12576 SCT2	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>		<input type="checkbox"/>
D17.23583 donor	<input type="radio"/>	<input type="radio"/>	<input checked="" type="radio"/>	<input type="radio"/>	<input type="radio"/>		<input type="checkbox"/>
D17.23584 host	<input type="radio"/>	<input checked="" type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>		<input type="checkbox"/>
D17.23585 SCT1	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>		<input type="checkbox"/>
D17.23585 SCT1 rep	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>		<input type="checkbox"/>

ASK ABOUT OUR OTHER APPLICATIONS

The MassARRAY is a versatile genetic analysis tool and not limited to a single application. Ask for information regarding our other applications including Oncology, PGx, hereditary genetic testing and sample integrity.

ORDERING INFORMATION

Catalog No.	Item	Chip format	# of Samples
13245D	Chimeric ID Panel Set - CPM (2x384)	384 CPM	96
13245	Chimeric ID Panel Set (2x384)	384	96
13169D	Chimeric ID Panel Set - CPM (10x384)	384 CPM	480
13169	Chimeric ID Panel Set (10x384)	384	480
13168F	Chimeric ID Panel Set - CPM (10x96)	96 CPM	120
13168	Chimeric ID Panel Set (10x96)	96	120

References

1. Pasquini MC, Wang Z: Current use and outcome of hematopoietic stem cell transplantation: CIBMTR Summary Slides, 2012. Available at: <http://www.cibmtr.org>.
2. NS Majhail, LW Mau, T Payton, E Denzen. National survey of blood and marrow transplant center personnel, infrastructure and models of care delivery. 2015. Available at: www.cibmtr.org.
3. JR Passweg et al. Hematopoietic stem cell transplantation in Europe 2014: more than 40 000 transplants annually. Bone Marrow Transplantation (2016) 51, 786-792.
4. Majhail et al. Prevalence of Hematopoietic Cell Transplant Survivors in the United States. Biol Blood Marror Transplant. 2013. 19(10): 1498-1501.

For Research Use Only. Not for use in diagnostic procedures.

Agena Bioscience, Inc.
4755 Eastgate Mall
San Diego, CA 92121
Phone: +1.858.882.2800

Order Desk: +1.858.202.9301
Order Desk Fax: +1.858.202.9220
orderdesk@agenabio.com
Web: agenabio.com

US +1.877.4.GENOME
EU +49.40.899676.0
AP +61.7.3088.1600
CN +86.21.6427.0566

